

For Research Use Only

Mouse VEGFA164 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00561-3

Capture Antibody Information

Catalog Number:
83572-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240415G1
Reactivity:
Mouse

Conjugate:
Unconjugated
Full name:
vascular endothelial growth factor A
Gene ID:
22339

Detection Antibody Information

Catalog Number:
83572-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240415A6
Reactivity:
Mouse
GenBank:
NM_001287057

Conjugate:
Unconjugated
Full name:
vascular endothelial growth factor A
Gene ID:
22339

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00561-3 targets VEGFA164 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse VEGFA164 Recombinant antibody, PBS Only (Capture) 83572-3-PBS (240415G1). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse VEGFA164 Recombinant antibody, PBS Only (Detector) 83572-1-PBS (240415A6). 100 µg. Concentration 1 mg/ml.

Alternative VEGFA164 matched antibody pairs: MP00561-1, MP00561-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00561-3, MOUSE VEGFA164 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83572-3-PBS. Detection antibody: 83572-1-PBS. Standard: Eg0277. Range: 0.625-80 ng/mL.