

For Research Use Only

BARHL2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00582-3

Capture Antibody Information

Catalog Number:
83596-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240558D3

Reactivity:
Human

Immunogen Catalog Number:
Ag21156

Conjugate:
Unconjugated

Full name:
BarH-like homeobox 2

Gene ID:
343472

Detection Antibody Information

Catalog Number:
83596-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240558F12

Reactivity:
Human

GenBank:
BC126439

Immunogen Catalog Number:
Ag21156

Conjugate:
Unconjugated

Full name:
BarH-like homeobox 2

Gene ID:
343472

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00582-3 targets BARHL2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BARHL2 Recombinant antibody, PBS Only (Capture) 83596-3-PBS (240558D3). 100 µg. Concentration 1 mg/ml.

Detection antibody: BARHL2 Recombinant antibody, PBS Only (Detector) 83596-2-PBS (240558F12). 100 µg. Concentration 1 mg/ml.

Alternative BARHL2 matched antibody pairs: MP00582-1, MP00582-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00582-3, BARHL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83596-3-PBS. Detection antibody: 83596-2-PBS. Standard: Ag21156. Range: 0.625-40 ng/mL