

For Research Use Only



# TADA1L Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00624-3

## Capture Antibody Information

Catalog Number:  
83641-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240637F8

Reactivity:  
Human

Immunogen Catalog Number:  
Ag14118

Conjugate:  
Unconjugated

Full name:  
transcriptional adaptor 1 (HFI1 homolog, yeast)-like

Gene ID:  
117143

## Detection Antibody Information

Catalog Number:  
83641-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240637H7

Reactivity:  
Human

GenBank:  
BC015401

Immunogen Catalog Number:  
Ag14118

Conjugate:  
Unconjugated

Full name:  
transcriptional adaptor 1 (HFI1 homolog, yeast)-like

Gene ID:  
117143

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00624-3 targets TADA1L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TADA1L Recombinant antibody, PBS Only (Capture) 83641-3-PBS (240637F8). 100 µg. Concentration 1 mg/ml.

Detection antibody: TADA1L Recombinant antibody, PBS Only (Detector) 83641-2-PBS (240637H7). 100 µg. Concentration 1 mg/ml.

Alternative TADA1L matched antibody pairs: MP00624-1, MP00624-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00624-3, TADA1L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83641-3-PBS. Detection antibody: 83641-2-PBS. Standard: Ag14118. Range: 0.625-80 ng/mL.