

For Research Use Only

PLA2G4E Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00638-3

Capture Antibody Information

Catalog Number:

83645-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240620G8

Reactivity:

Human

Immunogen Catalog Number:

Ag12808

Conjugate:

Unconjugated

Full name:

phospholipase A2, group IVE

Gene ID:

123745

Detection Antibody Information

Catalog Number:

83645-5-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240620E3

Reactivity:

Human

GenBank:

BC101584

Immunogen Catalog Number:

Ag12808

Conjugate:

Unconjugated

Full name:

phospholipase A2, group IVE

Gene ID:

123745

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00638-3 targets PLA2G4E in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PLA2G4E Recombinant antibody, PBS Only (Capture) 83645-2-PBS (240620G8). 100 µg. Concentration 1 mg/ml.

Detection antibody: PLA2G4E Recombinant antibody, PBS Only (Detector) 83645-5-PBS (240620E3). 100 µg. Concentration 1 mg/ml.

Alternative PLA2G4E matched antibody pairs: MP00638-1, MP00638-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00638-3, PLA2G4E Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83645-2-PBS. Detection antibody: 83645-5-PBS. Standard: Ag12808. Range: 0.625-40 ng/mL.