

For Research Use Only

KCNH7 Recombinant Matched Antibody Pair, PBS Only



Catalog Number: MP00649-3

Capture Antibody Information

Catalog Number:
83674-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240644B5

Reactivity:
Human

Immunogen Catalog Number:
Ag4532

Conjugate:
Unconjugated

Full name:
potassium voltage-gated channel,
subfamily H (eag-related), member 7

Gene ID:
90134

Detection Antibody Information

Catalog Number:
83674-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240644A11

Reactivity:
Human

GenBank:
BC035815

Immunogen Catalog Number:
Ag4532

Conjugate:
Unconjugated

Full name:
potassium voltage-gated channel,
subfamily H (eag-related), member 7

Gene ID:
90134

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-80 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP00649-3 targets KCNH7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KCNH7 Recombinant antibody, PBS Only (Capture/Detector) 83674-1-PBS (240644B5). 100 µg. Concentration 1 mg/ml.

Detection antibody: KCNH7 Recombinant antibody, PBS Only (Detector) 83674-3-PBS (240644A11). 100 µg. Concentration 1 mg/ml.

Alternative KCNH7 matched antibody pairs: MP00649-1, MP00649-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00649-3, KCNH7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83674-1-PBS. Detection antibody: 83674-3-PBS. Standard: Ag4532. Range: 1.25-80 ng/mL.