

For Research Use Only

Mouse CD155/PVR Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00705-1

Capture Antibody Information

Catalog Number:

83724-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240731B4

Reactivity:

mouse

Conjugate:

Unconjugated

Full name:

poliovirus receptor

Gene ID:

52118

Detection Antibody Information

Catalog Number:

83724-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240731A2

Reactivity:

mouse

GenBank:

NM_027514.2

Conjugate:

Unconjugated

Full name:

poliovirus receptor

Gene ID:

52118

Applications

Tested Applications:

Cytometric bead array

Range:

0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00705-1 targets CD155/PVR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD155/PVR Recombinant antibody, PBS Only (Capture) 83724-4-PBS (240731B4). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse CD155/PVR Recombinant antibody, PBS Only (Detector) 83724-1-PBS (240731A2). 100 µg. Concentration 1 mg/ml.

Alternative CD155/PVR matched antibody pairs: MP00705-2, MP00705-3, MP01053-1, MP01053-2, MP01053-3, MP50402-1, MP50402-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00705-1, MOUSE CD155/PVR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83724-4-PBS. Detection antibody: 83724-1-PBS. Standard: Eg0903. Range: 0.313-40 ng/mL