

# MIB1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00718-2

## Capture Antibody Information

**Catalog Number:**

83718-2-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

240676E6

**Reactivity:**

Human

**Immunogen Catalog Number:**

Ag2507

**Conjugate:**

Unconjugated

**Full name:**

mindbomb homolog 1 (Drosophila)

**Gene ID:**

57534

## Detection Antibody Information

**Catalog Number:**

83718-4-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

240676E10

**Reactivity:**

Human

**GenBank:**

BC022403

**Immunogen Catalog Number:**

Ag2507

**Conjugate:**

Unconjugated

**Full name:**

mindbomb homolog 1 (Drosophila)

**Gene ID:**

57534

## Applications

**Tested Applications:**

Cytometric bead array

**Range:**

1.25-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00718-2 targets MIB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MIB1 Recombinant antibody, PBS Only (Capture) 83718-2-PBS (240676E6). 100 µg. Concentration 1 mg/mL.

Detection antibody: MIB1 Recombinant antibody, PBS Only (Capture/Detector) 83718-4-PBS (240676E10). 100 µg. Concentration 1 mg/mL.

Alternative MIB1 matched antibody pairs: MP00718-1, MP00718-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology.

Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage buffer:**

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00718-2, MIB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83718-2-PBS. Detection antibody: 83718-4-PBS. Standard: Ag2507. Range: 1.25-80 ng/mL.