

For Research Use Only



Mouse Il6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00756-1

Capture Antibody Information

Catalog Number:
83795-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240701E8
Reactivity:
Mouse

Conjugate:
Unconjugated
Full name:
interleukin 6
Gene ID:
16193

Detection Antibody Information

Catalog Number:
83795-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240701B5
Reactivity:
Mouse
GenBank:
NM_031168

Conjugate:
Unconjugated
Full name:
interleukin 6
Gene ID:
16193

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00756-1 targets Il6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Il6 Recombinant antibody, PBS Only (Capture) 83795-1-PBS (240701E8). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Il6 Recombinant antibody, PBS Only (Detector) 83795-3-PBS (240701B5). 100 µg. Concentration 1 mg/mL.

Alternative Il6 matched antibody pairs: MP00687-1, MP00687-2, MP00687-3, MP00756-2, MP00756-3, MP50032-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00756-1, MOUSE IL-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83795-1-PBS. Detection antibody: 83795-3-PBS. Standard: Eg0973. Range: 1.25-40 ng/mL