

For Research Use Only

Islet 1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00771-4

Capture Antibody Information

Catalog Number:
83784-7-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240833C9

Reactivity:
human

Immunogen Catalog Number:
Ag33541

Conjugate:
Unconjugated

Full name:
ISL LIM homeobox 1

Gene ID:
3670

Detection Antibody Information

Catalog Number:
83784-6-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240833B2

Reactivity:
human

GenBank:
BC031213

Immunogen Catalog Number:
Ag33541

Conjugate:
Unconjugated

Full name:
ISL LIM homeobox 1

Gene ID:
3670

Applications

Tested Applications:
Sandwich ELISA

Range:
0.156-10 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00771-4 targets Islet 1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: Islet 1 Recombinant antibody, PBS Only (Capture) 83784-7-PBS (240833C9). 100 µg. Concentration 1 mg/ml.

Detection antibody: Islet 1 Recombinant antibody, PBS Only (Detector) 83784-6-PBS (240833B2). 100 µg. Concentration 1 mg/ml.

Alternative Islet 1 matched antibody pairs: MP00771-1, MP00771-2, MP00771-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00771-4, Human Islet 1 Recombinant Matched Antibody Pair - PBS only. 83784-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag33541. 83784-6-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.



The mean Islet 1 concentration was determined to be 1.05 ng/mL in SH-SY5Y cell extract based on a 3.40 mg/mL extract load.