

For Research Use Only

# DOT1L Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00836-2

## Capture Antibody Information

Catalog Number:  
83904-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240970F8

Reactivity:  
human

Immunogen Catalog Number:  
Ag27160

Conjugate:  
Unconjugated

Full name:  
DOT1-like, histone H3  
methyltransferase (*S. cerevisiae*)

Gene ID:  
84444

## Detection Antibody Information

Catalog Number:  
83904-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240970A5

Reactivity:  
human

GenBank:  
NM\_032482

Immunogen Catalog Number:  
Ag27160

Conjugate:  
Unconjugated

Full name:  
DOT1-like, histone H3  
methyltransferase (*S. cerevisiae*)

Gene ID:  
84444

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.56-100 ng/mL (Cytometric Bead  
Array)

Recommended Dilutions:  
It is recommended that this reagent  
should be titrated in each testing  
system to obtain optimal results.

## Product Information

MP00836-2 targets DOT1L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DOT1L Recombinant antibody, PBS Only (Capture) 83904-4-PBS (240970F8). 100 µg. Concentration 1 mg/ml.

Detection antibody: DOT1L Recombinant antibody, PBS Only (Detector) 83904-1-PBS (240970A5). 100 µg. Concentration 1 mg/ml.

Alternative DOT1L matched antibody pairs: MP00836-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00836-2, DOT1L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83904-4-PBS. Detection antibody: 83904-1-PBS. Standard: Ag27160. Range: 1.56-100 ng/mL