

IRAK3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00849-3

Capture Antibody Information

Catalog Number:

83876-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240975A5

Reactivity:

human

Immunogen Catalog Number:

Ag4959

Conjugate:

Unconjugated

Full name:

interleukin-1 receptor-associated kinase 3

Gene ID:

11213

Detection Antibody Information

Catalog Number:

83876-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240975C5

Reactivity:

human

GenBank:

BC057800

Immunogen Catalog Number:

Ag4959

Conjugate:

Unconjugated

Full name:

interleukin-1 receptor-associated kinase 3

Gene ID:

11213

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00849-3 targets IRAK3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IRAK3 Recombinant antibody, PBS Only (Capture) 83876-1-PBS (240975A5). 100 µg. Concentration 1 mg/mL.

Detection antibody: IRAK3 Recombinant antibody, PBS Only (Detector) 83876-2-PBS (240975C5). 100 µg. Concentration 1 mg/mL.

Alternative IRAK3 matched antibody pairs: MP00849-1, MP00849-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00849-3, IRAK3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83876-1-PBS. Detection antibody: 83876-2-PBS. Standard: Ag4959. Range: 0.625-80 ng/mL