

# GDF5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00881-2

## Capture Antibody Information

**Catalog Number:**

83988-2-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

240944E2

**Reactivity:**

human

**Immunogen Catalog Number:**

Ag6238

**Conjugate:**

Unconjugated

**Full name:**

growth differentiation factor 5

**Gene ID:**

8200

## Detection Antibody Information

**Catalog Number:**

83988-3-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

240944G6

**Reactivity:**

human

**GenBank:**

BC032495

**Immunogen Catalog Number:**

Ag6238

**Conjugate:**

Unconjugated

**Full name:**

growth differentiation factor 5

**Gene ID:**

8200

## Applications

**Tested Applications:**

Cytometric bead array

**Range:**

0.469-60 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00881-2 targets GDF5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GDF5 Recombinant antibody, PBS Only (Capture) 83988-2-PBS (240944E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: GDF5 Recombinant antibody, PBS Only (Capture/Detector) 83988-3-PBS (240944G6). 100 µg. Concentration 1 mg/mL.

Alternative GDF5 matched antibody pairs: MP00881-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology.

Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage buffer:**

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00881-2, GDF5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83988-2-PBS. Detection antibody: 83988-3-PBS. Standard: Ag6238. Range: 0.469-60 ng/mL