

For Research Use Only

CCDC22 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00890-2

Capture Antibody Information

Catalog Number:

83992-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241146G10

Reactivity:

human

Immunogen Catalog Number:

Ag9980

Conjugate:

Unconjugated

Full name:

coiled-coil domain containing 22

Gene ID:

28952

Detection Antibody Information

Catalog Number:

83992-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241146D11

Reactivity:

human

GenBank:

BC000972

Immunogen Catalog Number:

Ag9980

Conjugate:

Unconjugated

Full name:

coiled-coil domain containing 22

Gene ID:

28952

Applications

Tested Applications:

Cytometric bead array

Range:

0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00890-2 targets CCDC22 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCDC22 Recombinant antibody, PBS Only (Capture) 83992-3-PBS (241146G10). 100 µg. Concentration 1 mg/ml.

Detection antibody: CCDC22 Recombinant antibody, PBS Only (Detector) 83992-2-PBS (241146D11). 100 µg. Concentration 1 mg/ml.

Alternative CCDC22 matched antibody pairs: MP00890-1, MP00890-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology.

Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00890-2, CCDC22 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83992-3-PBS. Detection antibody: 83992-2-PBS. Standard: Ag9980. Range: 0.313-40 ng/mL