

For Research Use Only

QSOX1/Quiescin Q6 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00903-1

Capture Antibody Information

Catalog Number:

83952-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240946D5

Reactivity:

human

Immunogen Catalog Number:

Ag3396

Conjugate:

Unconjugated

Full name:

quiescin Q6 sulfhydryl oxidase 1

Gene ID:

5768

Detection Antibody Information

Catalog Number:

83952-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240946B3

Reactivity:

human

GenBank:

BC017692

Immunogen Catalog Number:

Ag3396

Conjugate:

Unconjugated

Full name:

quiescin Q6 sulfhydryl oxidase 1

Gene ID:

5768

Applications

Tested Applications:

Cytometric bead array

Range:

1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00903-1 targets QSOX1/Quiescin Q6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: QSOX1/Quiescin Q6 Recombinant antibody, PBS Only (Capture) 83952-2-PBS (240946D5). 100 µg. Concentration 1 mg/ml.

Detection antibody: QSOX1/Quiescin Q6 Recombinant antibody, PBS Only (Detector) 83952-1-PBS (240946B3). 100 µg. Concentration 1 mg/ml.

Alternative QSOX1/Quiescin Q6 matched antibody pairs: MP00903-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00903-1, QSOX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83952-2-PBS. Detection antibody: 83952-1-PBS. Standard: Ag3396. Range: 1.563-100 ng/mL.