

For Research Use Only

TLE3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00909-2

Capture Antibody Information

Catalog Number:
83955-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240680D11
Reactivity:
human
Immunogen Catalog Number:
Ag1930

Conjugate:
Unconjugated
Full name:
transducin-like enhancer of split 3 (E(sp1) homolog, Drosophila)
Gene ID:
7090

Detection Antibody Information

Catalog Number:
83955-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240680A9
Reactivity:
human
GenBank:
BC015729
Immunogen Catalog Number:
Ag1930

Conjugate:
Unconjugated
Full name:
transducin-like enhancer of split 3 (E(sp1) homolog, Drosophila)
Gene ID:
7090

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00909-2 targets TLE3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TLE3 Recombinant antibody, PBS Only (Capture/Detector) 83955-2-PBS (240680D11). 100 µg. Concentration 1 mg/ml.

Detection antibody: TLE3 Recombinant antibody, PBS Only (Detector) 83955-1-PBS (240680A9). 100 µg. Concentration 1 mg/ml.

Alternative TLE3 matched antibody pairs: MP00909-1, MP00909-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00909-2, TLE3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83955-2-PBS. Detection antibody: 83955-1-PBS. Standard: Ag1930. Range: 0.313-40 ng/mL