

For Research Use Only

# NFKB2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00950-3

## Capture Antibody Information

**Catalog Number:**  
84022-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241126D5  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag0662

**Conjugate:**  
Unconjugated  
**Full name:**  
nuclear factor of kappa light polypeptide gene enhancer in B-cells 2 (p49/p100)  
**Gene ID:**  
4791

## Detection Antibody Information

**Catalog Number:**  
84022-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241126F8  
**Reactivity:**  
human  
**GenBank:**  
BC002844  
**Immunogen Catalog Number:**  
Ag0662

**Conjugate:**  
Unconjugated  
**Full name:**  
nuclear factor of kappa light polypeptide gene enhancer in B-cells 2 (p49/p100)  
**Gene ID:**  
4791

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00950-3 targets NFKB2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NFKB2 Recombinant antibody, PBS Only (Capture) 84022-1-PBS (241126D5). 100 µg. Concentration 1 mg/ml.

Detection antibody: NFKB2 Recombinant antibody, PBS Only (Capture/Detector) 84022-4-PBS (241126F8). 100 µg. Concentration 1 mg/ml.

Alternative NFKB2 matched antibody pairs: MP00950-1, MP00950-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00950-3, NFkB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84022-1-PBS. Detection antibody: 84022-4-PBS. Standard: Ag0662. Range: Ag0662.