

For Research Use Only

# PCF11 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00996-1

## Capture Antibody Information

**Catalog Number:**  
84081-2-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
241213C11

**Reactivity:**  
human

**Immunogen Catalog Number:**  
Ag20263

**Conjugate:**  
Unconjugated

**Full name:**  
PCF11, cleavage and polyadenylation factor subunit, homolog (S. cerevisiae)

**Gene ID:**  
51585

## Detection Antibody Information

**Catalog Number:**  
84081-3-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
241213E8

**Reactivity:**  
human

**GenBank:**  
BC146778

**Immunogen Catalog Number:**  
Ag20263

**Conjugate:**  
Unconjugated

**Full name:**  
PCF11, cleavage and polyadenylation factor subunit, homolog (S. cerevisiae)

**Gene ID:**  
51585

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00996-1 targets PCF11 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PCF11 Recombinant antibody, PBS Only (Capture) 84081-2-PBS (241213C11). 100 µg. Concentration 1 mg/mL.

Detection antibody: PCF11 Recombinant antibody, PBS Only (Detector) 84081-3-PBS (241213E8). 100 µg. Concentration 1 mg/mL.

Alternative PCF11 matched antibody pairs: MP00996-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00996-1, PCF11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84081-2-PBS. Detection antibody: 84081-3-PBS. Standard: Ag20263. Range: 0.625-80 ng/mL