

PDLIM4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01021-2

Capture Antibody Information

Catalog Number:

84117-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241401E6

Reactivity:

human

Immunogen Catalog Number:

Ag13424

Conjugate:

Unconjugated

Full name:

PDZ and LIM domain 4

Gene ID:

8572

Detection Antibody Information

Catalog Number:

84117-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241401B3

Reactivity:

human

GenBank:

BC016765

Immunogen Catalog Number:

Ag13424

Conjugate:

Unconjugated

Full name:

PDZ and LIM domain 4

Gene ID:

8572

Applications

Tested Applications:

Cytometric bead array

Range:

0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01021-2 targets PDLIM4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PDLIM4 Recombinant antibody, PBS Only (Capture) 84117-2-PBS (241401E6). 100 µg. Concentration 1 mg/ml.

Detection antibody: PDLIM4 Recombinant antibody, PBS Only (Detector) 84117-1-PBS (241401B3). 100 µg. Concentration 1 mg/ml.

Alternative PDLIM4 matched antibody pairs: MP01021-1, MP01021-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MPO1021-2, PDLIM4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84117-2-PBS. Detection antibody: 84117-1-PBS. Standard: Ag13424. Range: 0.781-100 ng/mL.