

For Research Use Only



HESX1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01027-1

Capture Antibody Information

Catalog Number: 84125-4-PBS	Clone ID: 241430G8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: HESX homeobox 1
Isotype: IgG	Immunogen Catalog Number: Ag12348	Gene ID: 8820
Purification Method: Protein A purification		

Detection Antibody Information

Catalog Number: 84125-3-PBS	Clone ID: 241430G1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: HESX homeobox 1
Isotype: IgG	GenBank: BC093979	Gene ID: 8820
Purification Method: Protein A purification	Immunogen Catalog Number: Ag12348	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01027-1 targets HESX1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HESX1 Recombinant antibody, PBS Only (Capture/Detector) 84125-4-PBS (241430G8). 100 µg. Concentration 1 mg/ml.

Detection antibody: HESX1 Recombinant antibody, PBS Only (Detector) 84125-3-PBS (241430G1). 100 µg. Concentration 1 mg/ml.

Alternative HESX1 matched antibody pairs: MP01027-2, MP01027-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01027-1, HESX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84125-4-PBS. Detection antibody: 84125-3-PBS. Standard: Ag12348. Range: 0.781-100 ng/mL.