

For Research Use Only

# Mouse HER2/ErbB2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01039-1

## Capture Antibody Information

**Catalog Number:**  
84128-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241366G10  
**Reactivity:**  
mouse

**Conjugate:**  
Unconjugated  
**Full name:**  
v-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian)  
**Gene ID:**  
13866

## Detection Antibody Information

**Catalog Number:**  
84128-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241366D1  
**Reactivity:**  
mouse  
**GenBank:**  
NM\_001003817.1

**Conjugate:**  
Unconjugated  
**Full name:**  
v-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian)  
**Gene ID:**  
13866

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.469-60 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01039-1 targets HER2/ErbB2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: Mouse HER2/ErbB2 Recombinant antibody, PBS Only (Capture) 84128-3-PBS (241366G10). 100 µg. Concentration 1 mg/ml.  
Detection antibody: Mouse HER2/ErbB2 Recombinant antibody, PBS Only (Detector) 84128-2-PBS (241366D1). 100 µg. Concentration 1 mg/ml.  
Alternative HER2/ErbB2 matched antibody pairs: MP01039-2, MP01039-3, MP01039-4, MP01080-1, MP01080-2, MP01080-3  
Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.  
Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01039-1, MOUSE Her2/Erbb2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84128-3-PBS. Detection antibody: 84128-2-PBS. Standard: Eg0803. Range: 0.469-60 ng/mL