

For Research Use Only

Mouse HER2/ErbB2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01039-2

Capture Antibody Information

Catalog Number:
84128-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241366C3
Reactivity:
mouse

Conjugate:
Unconjugated
Full name:
v-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian)
Gene ID:
13866

Detection Antibody Information

Catalog Number:
84128-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241366D1
Reactivity:
mouse
GenBank:
NM_001003817.1

Conjugate:
Unconjugated
Full name:
v-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian)
Gene ID:
13866

Applications

Tested Applications:
Cytometric bead array

Range:
0.469-60 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01039-2 targets HER2/ErbB2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: Mouse HER2/ErbB2 Recombinant antibody, PBS Only (Capture) 84128-1-PBS (241366C3). 100 µg. Concentration 1 mg/mL.
Detection antibody: Mouse HER2/ErbB2 Recombinant antibody, PBS Only (Detector) 84128-2-PBS (241366D1). 100 µg. Concentration 1 mg/mL.
Alternative HER2/ErbB2 matched antibody pairs: MP01039-1, MP01039-3, MP01039-4, MP01080-1, MP01080-2, MP01080-3
Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01039-2, MOUSE Her2/ErbB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84128-1-PBS. Detection antibody: 84128-2-PBS. Standard: Eg0803. Range: 0.469-60 ng/mL