

For Research Use Only

Rat gp130/IL6ST Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01049-2

Capture Antibody Information

Catalog Number:

84113-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241340G7

Reactivity:

rat

Conjugate:

Unconjugated

Full name:

interleukin 6 signal transducer

Gene ID:

25205

Detection Antibody Information

Catalog Number:

84113-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241340C2

Reactivity:

rat

GenBank:

NM_001008725.3

Conjugate:

Unconjugated

Full name:

interleukin 6 signal transducer

Gene ID:

25205

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01049-2 targets gp130/IL6ST in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Rat gp130/IL6ST Recombinant antibody, PBS Only (Capture) 84113-4-PBS (241340G7). 100 µg. Concentration 1 mg/ml.

Detection antibody: Rat gp130/IL6ST Recombinant antibody, PBS Only (Detector) 84113-1-PBS (241340C2). 100 µg. Concentration 1 mg/ml.

Alternative gp130/IL6ST matched antibody pairs: MP00324-1, MP00324-2, MP01049-1, MP01049-3, MP50086-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology.

Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01049-2, RAT gp130 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84113-4-PBS. Detection antibody: 84113-1-PBS. Standard: Eg0857. Range: 0.625-80 ng/mL.