

For Research Use Only

# Mouse CD27/TNFRSF7 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01069-3

## Capture Antibody Information

Catalog Number:  
84157-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241378A3

Reactivity:  
mouse

Conjugate:  
Unconjugated

Full name:  
CD27 antigen

Gene ID:  
21940

## Detection Antibody Information

Catalog Number:  
84157-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241378E11

Reactivity:  
mouse

GenBank:  
NM\_001033126.2

Conjugate:  
Unconjugated

Full name:  
CD27 antigen

Gene ID:  
21940

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01069-3 targets CD27/TNFRSF7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD27/TNFRSF7 Recombinant antibody, PBS Only (Capture) 84157-1-PBS (241378A3). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse CD27/TNFRSF7 Recombinant antibody, PBS Only (Capture/Detector) 84157-4-PBS (241378E11). 100 µg. Concentration 1 mg/ml.

Alternative CD27/TNFRSF7 matched antibody pairs: MP01069-1, MP01069-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01069-3, MOUSE CD27/TNFRSF7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84157-1-PBS. Detection antibody: 84157-4-PBS. Standard: Eg1311. Range: 0.781-50 ng/mL