

For Research Use Only

LIMP2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01086-4

Capture Antibody Information

Catalog Number:
84193-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241468F9

Reactivity:
human

Immunogen Catalog Number:
Ag25908

Conjugate:
Unconjugated

Full name:
scavenger receptor class B, member 2

Gene ID:
950

Detection Antibody Information

Catalog Number:
84193-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241468G1

Reactivity:
human

GenBank:
BC021892

Immunogen Catalog Number:
Ag25908

Conjugate:
Unconjugated

Full name:
scavenger receptor class B, member 2

Gene ID:
950

Applications

Tested Applications:
Sandwich ELISA

Range:
3.91-250 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01086-4 targets LIMP2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: LIMP2 Recombinant antibody, PBS Only (Capture) 84193-4-PBS (241468F9). 100 µg. Concentration 1 mg/mL.

Detection antibody: LIMP2 Recombinant antibody, PBS Only (Detector) 84193-3-PBS (241468G1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01086-4, Human LIMP2 Recombinant Matched Antibody Pair - PBS only. 84193-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag25908. 84193-3-PBS was HRP conjugated as the detection antibody. Range: 3.91-250 ng/mL.



The mean LIMP2 concentration was determined to be 271.42 ng/mL in mouse liver cell extract based on a 3.70 mg/mL extract load.