

For Research Use Only

4-1BBL/TNFSF9 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01091-1

Capture Antibody Information

Catalog Number:
84185-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241386B11

Reactivity:
human

Conjugate:
Unconjugated

Full name:
tumor necrosis factor (ligand)
superfamily, member 9

Gene ID:
8744

Detection Antibody Information

Catalog Number:
84185-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241386D5

Reactivity:
human

GenBank:
BC104805

Conjugate:
Unconjugated

Full name:
tumor necrosis factor (ligand)
superfamily, member 9

Gene ID:
8744

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01091-1 targets 4-1BBL/TNFSF9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: 4-1BBL/TNFSF9 Recombinant antibody, PBS Only (Capture) 84185-2-PBS (241386B11). 100 µg. Concentration 1 mg/ml.

Detection antibody: 4-1BBL/TNFSF9 Recombinant antibody, PBS Only (Detector) 84185-3-PBS (241386D5). 100 µg. Concentration 1 mg/ml.

Alternative 4-1BBL/TNFSF9 matched antibody pairs: MP01091-2, MP50011-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01091-1, 4-1BBL/TNFSF9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84185-2-PBS. Detection antibody: 84185-3-PBS. Standard: Eg0010. Range: 0.156-20 ng/mL.