

For Research Use Only

PDCD4 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01098-1

Capture Antibody Information

Catalog Number:
84162-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241463E3

Reactivity:
human

Immunogen Catalog Number:
Ag3288

Conjugate:
Unconjugated

Full name:
programmed cell death 4 (neoplastic transformation inhibitor)

Gene ID:
27250

Detection Antibody Information

Catalog Number:
84162-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241463A8

Reactivity:
human

GenBank:
BC026104

Immunogen Catalog Number:
Ag3288

Conjugate:
Unconjugated

Full name:
programmed cell death 4 (neoplastic transformation inhibitor)

Gene ID:
27250

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01098-1 targets PDCD4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PDCD4 Recombinant antibody, PBS Only (Capture) 84162-3-PBS (241463E3). 100 µg. Concentration 1 mg/ml.

Detection antibody: PDCD4 Recombinant antibody, PBS Only (Capture/Detector) 84162-1-PBS (241463A8). 100 µg. Concentration 1 mg/ml.

Alternative PDCD4 matched antibody pairs: MP01098-2, MP01098-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01098-1, PDCD4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84162-3-PBS. Detection antibody: 84162-1-PBS. Standard: Ag3288. Range: 0.781-100 ng/mL.