

For Research Use Only

C19orf33 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01103-2

Capture Antibody Information

Catalog Number:
84148-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241220D10

Reactivity:
human

Immunogen Catalog Number:
Ag22218

Conjugate:
Unconjugated

Full name:
chromosome 19 open reading frame
33

Gene ID:
64073

Detection Antibody Information

Catalog Number:
84148-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241220C6

Reactivity:
human

GenBank:
BC060319

Immunogen Catalog Number:
Ag22218

Conjugate:
Unconjugated

Full name:
chromosome 19 open reading frame
33

Gene ID:
64073

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01103-2 targets C19orf33 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C19orf33 Recombinant antibody, PBS Only (Capture) 84148-4-PBS (241220D10). 100 µg. Concentration 1 mg/ml.

Detection antibody: C19orf33 Recombinant antibody, PBS Only (Detector) 84148-3-PBS (241220C6). 100 µg. Concentration 1 mg/ml.

Alternative C19orf33 matched antibody pairs: MP01103-1, MP01103-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01103-2, C19orf33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84148-4-PBS. Detection antibody: 84148-3-PBS. Standard: Ag22218. Range: 0.313-40 ng/mL.