

For Research Use Only

# YTHDF1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01146-4

## Capture Antibody Information

Catalog Number:  
80876-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241445F12

Reactivity:  
human

Immunogen Catalog Number:  
Ag11509

Conjugate:  
Unconjugated

Full name:  
YTH domain family, member 1

Gene ID:  
54915

## Detection Antibody Information

Catalog Number:  
80876-7-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241445A1

Reactivity:  
human

GenBank:  
BC050284

Immunogen Catalog Number:  
Ag11509

Conjugate:  
Unconjugated

Full name:  
YTH domain family, member 1

Gene ID:  
54915

## Applications

Tested Applications:  
Sandwich ELISA

Range:  
78.1-5000 pg/mL (Sandwich ELISA)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01146-4 targets YTHDF1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: YTHDF1 Recombinant antibody, PBS Only (Capture) 80876-4-PBS (241445F12). 100 µg. Concentration 1 mg/ml.

Detection antibody: Ag11509 Recombinant antibody, PBS Only (Detector) 80876-7-PBS (241445A1). 100 µg. Concentration 1 mg/ml.

Alternative YTHDF1 matched antibody pairs: MP01146-1, MP01146-2, MP01146-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP01146-4, Human YTHDF1 Recombinant Matched Antibody Pair - PBS only. 80876-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag11509. 80876-7-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean YTHDF1 concentration was determined to be 22.0 ng/mL in HEK-293 cell extract based on a 1.5 mg/mL extract load.