

For Research Use Only

Noggin Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01174-4

Capture Antibody Information

Catalog Number:
84283-6-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241511E5

Reactivity:
human

Immunogen Catalog Number:
Ag32755

Conjugate:
Unconjugated

Full name:
noggin

Gene ID:
9241

Detection Antibody Information

Catalog Number:
84283-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241511D9

Reactivity:
human

GenBank:
BC034027

Immunogen Catalog Number:
Ag32755

Conjugate:
Unconjugated

Full name:
noggin

Gene ID:
9241

Applications

Tested Applications:
Sandwich ELISA

Range:
0.195-12.5 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01174-4 targets Noggin in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: NOG Recombinant antibody, PBS Only (Capture) 84283-6-PBS (241511E5). 100 µg. Concentration 1 mg/ml.

Detection antibody: NOG Recombinant antibody, PBS Only (Detector) 84283-3-PBS (241511D9). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01174-4, Human Noggin Recombinant Matched Antibody Pair - PBS only. 84283-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag32755. 84283-3-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL.



The mean Noggin concentration was determined to be 79.96 ng/mL in mouse brain tissue extract based on a 2.0 mg/mL extract load.