

For Research Use Only

IVNS1ABP Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01182-2

Capture Antibody Information

Catalog Number:
84268-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241587H4

Reactivity:
human

Immunogen Catalog Number:
Ag6412

Conjugate:
Unconjugated

Full name:
influenza virus NS1A binding protein

Gene ID:
10625

Detection Antibody Information

Catalog Number:
84268-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241587B1

Reactivity:
human

GenBank:
BC067739

Immunogen Catalog Number:
Ag6412

Conjugate:
Unconjugated

Full name:
influenza virus NS1A binding protein

Gene ID:
10625

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01182-2 targets IVNS1ABP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IVNS1ABP Recombinant antibody, PBS Only (Capture) 84268-2-PBS (241587H4). 100 µg. Concentration 1 mg/ml.

Detection antibody: IVNS1ABP Recombinant antibody, PBS Only (Detector) 84268-1-PBS (241587B1). 100 µg. Concentration 1 mg/ml.

Alternative IVNS1ABP matched antibody pairs: MP01182-1, MP01182-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01182-2, IVNS1ABP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84268-2-PBS. Detection antibody: 84268-1-PBS. Standard: Ag6412. Range: 1.25-40 ng/mL.