

For Research Use Only

NME3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01200-1

Capture Antibody Information

Catalog Number:
84278-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241585H2
Reactivity:
human
Immunogen Catalog Number:
Ag7284

Conjugate:
Unconjugated
Full name:
non-metastatic cells 3, protein expressed in
Gene ID:
4832

Detection Antibody Information

Catalog Number:
84278-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241585D6
Reactivity:
human
GenBank:
BC000250
Immunogen Catalog Number:
Ag7284

Conjugate:
Unconjugated
Full name:
non-metastatic cells 3, protein expressed in
Gene ID:
4832

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01200-1 targets NME3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NME3 Recombinant antibody, PBS Only (Capture) 84278-2-PBS (241585H2). 100 µg. Concentration 1 mg/mL.

Detection antibody: NME3 Recombinant antibody, PBS Only (Detector) 84278-1-PBS (241585D6). 100 µg. Concentration 1 mg/mL.

Alternative NME3 matched antibody pairs: MP01200-2, MP01200-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01200-1, NME3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84278-2-PBS. Detection antibody: 84278-1-PBS. Standard: Ag7284. Range: 0.313-40 ng/mL