

For Research Use Only



Mouse TNFRSF9/CD137 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01216-2

Capture Antibody Information

Catalog Number:
84315-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241647D8
Reactivity:
mouse

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor
superfamily, member 9
Gene ID:
21942

Detection Antibody Information

Catalog Number:
84315-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241647G11
Reactivity:
mouse
GenBank:
NM_001077509.1

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor
superfamily, member 9
Gene ID:
21942

Applications

Tested Applications:
Cytometric bead array

Range:
0.078-10 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01216-2 targets TNFRSF9/CD137 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse 4-1BB/TNFRSF9 Recombinant antibody, PBS Only (Capture) 84315-1-PBS (241647D8). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse 4-1BB/TNFRSF9 Recombinant antibody, PBS Only (Detector) 84315-3-PBS (241647G11). 100 µg. Concentration 1 mg/ml.

Alternative TNFRSF9/CD137 matched antibody pairs: MP00870-1, MP00870-2, MP50110-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

**This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.**

Selected Validation Data



Cytometric bead array standard curve of MP01216-2, MOUSE 4-1BB/TNFRSF9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84315-1-PBS. Detection antibody: 84315-3-PBS. Standard: Eg1627. Range: 0.078-10 ng/mL.