

For Research Use Only

# POLR3C Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01234-1

## Capture Antibody Information

Catalog Number:  
84367-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241402E12

Reactivity:  
human

Immunogen Catalog Number:  
Ag27311

Conjugate:  
Unconjugated

Full name:  
polymerase (RNA) III (DNA directed)  
polypeptide C (62kD)

Gene ID:  
10623

## Detection Antibody Information

Catalog Number:  
84367-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241402C9

Reactivity:  
human

GenBank:  
BC002586

Immunogen Catalog Number:  
Ag27311

Conjugate:  
Unconjugated

Full name:  
polymerase (RNA) III (DNA directed)  
polypeptide C (62kD)

Gene ID:  
10623

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01234-1 targets POLR3C in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: POLR3C Recombinant antibody, PBS Only (Capture/Detector) 84367-4-PBS (241402E12). 100 µg. Concentration 1 mg/ml.

Detection antibody: POLR3C Recombinant antibody, PBS Only (Detector) 84367-3-PBS (241402C9). 100 µg. Concentration 1 mg/ml.

Alternative POLR3C matched antibody pairs: MP01234-2, MP01234-3, MP01234-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01234-1, POLR3C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84367-4-PBS. Detection antibody: 84367-3-PBS. Standard: Ag27311. Range: 0.781-100 ng/mL