

For Research Use Only

IKKA Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01243-1

Capture Antibody Information

Catalog Number:
84372-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241579G10
Reactivity:
human
Immunogen Catalog Number:
Ag26573

Conjugate:
Unconjugated
Full name:
conserved helix-loop-helix
ubiquitous kinase
Gene ID:
1147

Detection Antibody Information

Catalog Number:
84372-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241579G3
Reactivity:
human
GenBank:
NM_001278
Immunogen Catalog Number:
Ag26573

Conjugate:
Unconjugated
Full name:
conserved helix-loop-helix
ubiquitous kinase
Gene ID:
1147

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01243-1 targets IKKA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IKKA Recombinant antibody, PBS Only (Capture) 84372-3-PBS (241579G10). 100 µg. Concentration 1 mg/ml.

Detection antibody: IKKA Recombinant antibody, PBS Only (Detector) 84372-2-PBS (241579G3). 100 µg. Concentration 1 mg/ml.

Alternative IKKA matched antibody pairs: MP01243-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01243-1, IKKA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84372-3-PBS. Detection antibody: 84372-2-PBS. Standard: Ag26573. Range: 0.781-100 ng/mL.