

For Research Use Only

CYP2R1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01251-1

Capture Antibody Information

Catalog Number:
84378-4-RR
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241072F5
Reactivity:
human
Immunogen Catalog Number:
Ag12873

Conjugate:
Unconjugated
Full name:
cytochrome P450, family 2, subfamily R, polypeptide 1
Gene ID:
120227

Detection Antibody Information

Catalog Number:
84378-2-RR
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241072E2
Reactivity:
human
GenBank:
BC104909
Immunogen Catalog Number:
Ag12873

Conjugate:
Unconjugated
Full name:
cytochrome P450, family 2, subfamily R, polypeptide 1
Gene ID:
120227

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01251-1 targets CYP2R1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: CYP2R1 Recombinant antibody 84378-4-RR (241072F5). 100 µg. Concentration 1 mg/mL.
Detection antibody: CYP2R1 Recombinant antibody 84378-2-RR (241072E2). 100 µg. Concentration 1 mg/mL.
Alternative CYP2R1 matched antibody pairs: MP01251-2, MP01251-3
Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01251-1, CYP2R1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84377-2-PBS. Detection antibody: 84377-1-PBS. Standard: Ag12873. Range: 0.781-100 ng/mL.