

For Research Use Only

ACE Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01285-2

Capture Antibody Information

Catalog Number:
84407-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241763D8
Reactivity:
human

Conjugate:
Unconjugated
Full name:
angiotensin I converting enzyme (peptidyl-dipeptidase A) 1
Gene ID:
1636

Detection Antibody Information

Catalog Number:
84407-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241763C6
Reactivity:
human
GenBank:
NM_000789.4

Conjugate:
Unconjugated
Full name:
angiotensin I converting enzyme (peptidyl-dipeptidase A) 1
Gene ID:
1636

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01285-2 targets ACE in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ACE Recombinant antibody, PBS Only (Capture) 84407-3-PBS (241763D8). 100 µg. Concentration 1 mg/ml.

Detection antibody: ACE Recombinant antibody, PBS Only (Detector) 84407-2-PBS (241763C6). 100 µg. Concentration 1 mg/ml.

Alternative ACE matched antibody pairs: MP01285-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01285-2, ACE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84407-3-PBS. Detection antibody: 84407-2-PBS. Standard: Eg2298. Range: 1.563-100 ng/mL.