

For Research Use Only

Mouse CEACAM-1/CD66a Recombinant Matched Antibody Pair, PBS Only



Catalog Number: **MP01297-3**

Capture Antibody Information

Catalog Number:
84409-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241856E3
Reactivity:
mouse

Conjugate:
Unconjugated
Full name:
carcinoembryonic antigen-related
cell adhesion molecule 1
Gene ID:
26365

Detection Antibody Information

Catalog Number:
84409-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241856H5
Reactivity:
mouse
GeneBank:

Conjugate:
Unconjugated
Full name:
carcinoembryonic antigen-related
cell adhesion molecule 1
Gene ID:
26365

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01297-3 targets CEACAM-1/CD66a in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CEACAM-1/CD66a Recombinant antibody, PBS Only (Capture) 84409-2-PBS (241856E3). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse CEACAM-1/CD66a Recombinant antibody, PBS Only (Capture/Detector) 84409-4-PBS (241856H5). 100 µg. Concentration 1 mg/ml.

Alternative CEACAM-1/CD66a matched antibody pairs: MP00520-1, MP01297-1, MP01297-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

**This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.**

Selected Validation Data



Cytometric bead array standard curve of MP01297-3, MOUSE CEACAM1/CD66a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84409-2-PBS. Detection antibody: 84409-4-PBS. Standard: Eg1423. Range: 0.313-40 ng/mL.