

For Research Use Only

CCL20/MIP-3 alpha Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01301-2

Capture Antibody Information

Catalog Number:
84413-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241694A5

Reactivity:
human

Conjugate:
Unconjugated

Full name:
chemokine (C-C motif) ligand 20

Gene ID:
6364

Detection Antibody Information

Catalog Number:
84413-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241694B2

Reactivity:
human

GenBank:
NM_004591.3

Conjugate:
Unconjugated

Full name:
chemokine (C-C motif) ligand 20

Gene ID:
6364

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01301-2 targets CCL20/MIP-3 alpha in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCL20/MIP-3 alpha Recombinant antibody, PBS Only (Capture) 84413-1-PBS (241694A5). 100 µg. Concentration 1 mg/ml.

Detection antibody: CCL20/MIP-3 alpha Recombinant antibody, PBS Only (Detector) 84413-2-PBS (241694B2). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01301-2, CCL20/MIP-3 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84413-1-PBS. Detection antibody: 84413-2-PBS. Standard: Eg2398. Range: 0.781-100 ng/mL.