

For Research Use Only

# Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01329-3

## Capture Antibody Information

**Catalog Number:**  
84460-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241725C7  
**Reactivity:**  
human

**Conjugate:**  
Unconjugated  
**Full name:**  
poliovirus receptor-related 2  
(herpesvirus entry mediator B)  
**Gene ID:**  
5819

## Detection Antibody Information

**Catalog Number:**  
84460-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241725A10  
**Reactivity:**  
human  
**GenBank:**  
NM\_001042724.2

**Conjugate:**  
Unconjugated  
**Full name:**  
poliovirus receptor-related 2  
(herpesvirus entry mediator B)  
**Gene ID:**  
5819

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01329-3 targets Nectin-2/CD112 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Nectin-2/CD112 Recombinant antibody, PBS Only (Capture) 84460-3-PBS (241725C7). 100 µg. Concentration 1 mg/ml.

Detection antibody: Nectin-2/CD112 Recombinant antibody, PBS Only (Detector) 84460-1-PBS (241725A10). 100 µg. Concentration 1 mg/ml.

Alternative Nectin-2/CD112 matched antibody pairs: MP01292-1, MP01292-2, MP01292-3, MP01329-1, MP01329-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01329-3, Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84460-3-PBS. Detection antibody: 84460-1-PBS. Standard: Eg1670. Range: 0.781-100 ng/mL.