

For Research Use Only

TLR2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01332-3

Capture Antibody Information

Catalog Number:

84462-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241681D11

Reactivity:

human

Conjugate:

Unconjugated

Full name:

toll-like receptor 2

Gene ID:

7097

Detection Antibody Information

Catalog Number:

84462-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241681E10

Reactivity:

human

GenBank:

BC033756

Conjugate:

Unconjugated

Full name:

toll-like receptor 2

Gene ID:

7097

Applications

Tested Applications:

Cytometric bead array

Range:

0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01332-3 targets TLR2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TLR2 Recombinant antibody, PBS Only (Capture) 84462-2-PBS (241681D11). 100 µg. Concentration 1 mg/mL.

Detection antibody: TLR2 Recombinant antibody, PBS Only (Detector) 84462-3-PBS (241681E10). 100 µg. Concentration 1 mg/mL.

Alternative TLR2 matched antibody pairs: MP01332-1, MP01332-2, MP50388-1, MP50388-2, MP50388-3, MP50388-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology.

Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01332-3, TLR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84462-2-PBS. Detection antibody: 84462-3-PBS. Standard: RP00009. Range: 0.781-100 ng/mL.