

For Research Use Only

LMO7 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01336-2

Capture Antibody Information

Catalog Number:
84495-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241781D12

Reactivity:
human

Immunogen Catalog Number:
Ag30016

Conjugate:
Unconjugated

Full name:
LIM domain 7

Gene ID:
4008

Detection Antibody Information

Catalog Number:
84495-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241781C9

Reactivity:
human

GenBank:
NM_005358

Immunogen Catalog Number:
Ag30016

Conjugate:
Unconjugated

Full name:
LIM domain 7

Gene ID:
4008

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01336-2 targets LMO7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LMO7 Recombinant antibody, PBS Only (Capture) 84495-2-PBS (241781D12). 100 µg. Concentration 1 mg/ml.

Detection antibody: LMO7 Recombinant antibody, PBS Only (Capture/Detector) 84495-1-PBS (241781C9). 100 µg. Concentration 1 mg/ml.

Alternative LMO7 matched antibody pairs: MP01336-1, MP01336-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01336-2, LMO7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84495-2-PBS. Detection antibody: 84495-1-PBS. Standard: Ag30016. Range: 0.625-80 ng/mL