

For Research Use Only

# FKBP11 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01342-3

## Capture Antibody Information

Catalog Number:  
84478-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241912H7

Reactivity:  
human

Immunogen Catalog Number:  
Ag11368

Conjugate:  
Unconjugated

Full name:  
FK506 binding protein 11, 19 kDa

Gene ID:  
51303

## Detection Antibody Information

Catalog Number:  
84478-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241912B1

Reactivity:  
human

GenBank:  
BC027973

Immunogen Catalog Number:  
Ag11368

Conjugate:  
Unconjugated

Full name:  
FK506 binding protein 11, 19 kDa

Gene ID:  
51303

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01342-3 targets FKBP11 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FKBP11 Recombinant antibody, PBS Only (Capture/Detector) 84478-4-PBS (241912H7). 100 µg. Concentration 1 mg/ml.

Detection antibody: FKBP11 Recombinant antibody, PBS Only (Detector) 84478-1-PBS (241912B1). 100 µg. Concentration 1 mg/ml.

Alternative FKBP11 matched antibody pairs: MP01342-1, MP01342-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01342-3, FKBP11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84478-4-PBS. Detection antibody: 84478-1-PBS. Standard: Ag11368. Range: 0.313-40 ng/mL