

For Research Use Only

# LASS2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01385-2

## Capture Antibody Information

Catalog Number:  
84514-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241898C5

Reactivity:  
human

Immunogen Catalog Number:  
Ag14151

Conjugate:  
Unconjugated

Full name:  
LAG1 homolog, ceramide synthase 2

Gene ID:  
29956

## Detection Antibody Information

Catalog Number:  
84514-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241898D10

Reactivity:  
human

GenBank:  
BC010032

Immunogen Catalog Number:  
Ag14151

Conjugate:  
Unconjugated

Full name:  
LAG1 homolog, ceramide synthase 2

Gene ID:  
29956

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01385-2 targets LASS2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LASS2 Recombinant antibody, PBS Only (Capture/Detector) 84514-1-PBS (241898C5). 100 µg. Concentration 1 mg/ml.

Detection antibody: LASS2 Recombinant antibody, PBS Only (Detector) 84514-3-PBS (241898D10). 100 µg. Concentration 1 mg/ml.

Alternative LASS2 matched antibody pairs: MP01385-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01385-2, LASS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84514-1-PBS. Detection antibody: 84514-3-PBS. Standard: Ag14151. Range: 0.156-20 ng/mL