

For Research Use Only

FADD Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01426-1

Capture Antibody Information

Catalog Number:
84619-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241976E11
Reactivity:
human
Immunogen Catalog Number:
Ag6701

Conjugate:
Unconjugated
Full name:
Fas (TNFRSF6)-associated via death domain
Gene ID:
8772

Detection Antibody Information

Catalog Number:
84619-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241976E2
Reactivity:
human
GenBank:
BC000334
Immunogen Catalog Number:
Ag6701

Conjugate:
Unconjugated
Full name:
Fas (TNFRSF6)-associated via death domain
Gene ID:
8772

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01426-1 targets FADD in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FADD Recombinant antibody, PBS Only (Capture) 84619-4-PBS (241976E11). 100 µg. Concentration 1 mg/ml.

Detection antibody: FADD Recombinant antibody, PBS Only (Detector) 84619-3-PBS (241976E2). 100 µg. Concentration 1 mg/ml.

Alternative FADD matched antibody pairs: MP01426-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01426-1, FADD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84619-4-PBS. Detection antibody: 84619-3-PBS. Standard: Ag6701. Range: 0.781-100 ng/mL.