

For Research Use Only

# ATP8A2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01468-2

## Capture Antibody Information

Catalog Number:  
84653-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242105G7

Reactivity:  
human

Immunogen Catalog Number:  
Ag36276

Conjugate:  
Unconjugated

Full name:  
ATPase, aminophospholipid transporter-like, class I, type 8A, member 2

Gene ID:  
51761

## Detection Antibody Information

Catalog Number:  
84653-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242105H2

Reactivity:  
human

GenBank:  
NM\_001313741

Immunogen Catalog Number:  
Ag36276

Conjugate:  
Unconjugated

Full name:  
ATPase, aminophospholipid transporter-like, class I, type 8A, member 2

Gene ID:  
51761

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.25-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01468-2 targets ATP8A2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ATP8A2 Recombinant antibody, PBS Only (Capture) 84653-2-PBS (242105G7). 100 µg. Concentration 1 mg/ml.

Detection antibody: ATP8A2 Recombinant antibody, PBS Only (Detector) 84653-3-PBS (242105H2). 100 µg. Concentration 1 mg/ml.

Alternative ATP8A2 matched antibody pairs: MP01468-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01468-2, ATP8A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84653-2-PBS. Detection antibody: 84653-3-PBS. Standard: Ag36276. Range: 1.25-80 ng/mL.