

For Research Use Only

SEMA7A Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01498-1

Capture Antibody Information

Catalog Number:
84682-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241814B12
Reactivity:
human

Conjugate:
Unconjugated
Full name:
semaphorin 7A, GPI membrane anchor (John Milton Hagen blood group)
Gene ID:
8482

Detection Antibody Information

Catalog Number:
84682-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241814A12
Reactivity:
human
GenBank:
NM_003612.5

Conjugate:
Unconjugated
Full name:
semaphorin 7A, GPI membrane anchor (John Milton Hagen blood group)
Gene ID:
8482

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01498-1 targets SEMA7A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Semaphorin 7A Recombinant antibody, PBS Only (Capture) 84682-2-PBS (241814B12). 100 µg. Concentration 1 mg/ml.

Detection antibody: Semaphorin 7A Recombinant antibody, PBS Only (Detector) 84682-1-PBS (241814A12). 100 µg. Concentration 1 mg/ml.

Alternative SEMA7A matched antibody pairs: MP50321-1, MP50321-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01498-1, Semaphorin 7A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84682-2-PBS. Detection antibody: 84682-1-PBS. Standard: Eg2005. Range: 0.781-100 ng/mL.