

For Research Use Only



TMIGD2/CD28H Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01528-2

Capture Antibody Information

Catalog Number:
84727-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241807A2
Reactivity:
human

Conjugate:
Unconjugated
Full name:
transmembrane and immunoglobulin domain containing 2
Gene ID:
126259

Detection Antibody Information

Catalog Number:
84727-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241807C10
Reactivity:
human
GenBank:
BC015655

Conjugate:
Unconjugated
Full name:
transmembrane and immunoglobulin domain containing 2
Gene ID:
126259

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01528-2 targets TMIGD2/CD28H in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TMIGD2/CD28H Recombinant antibody, PBS Only (Capture) 84727-1-PBS (241807A2). 100 µg. Concentration 1 mg/mL.

Detection antibody: TMIGD2/CD28H Recombinant antibody, PBS Only (Detector) 84727-3-PBS (241807C10). 100 µg. Concentration 1 mg/mL.

Alternative TMIGD2/CD28H matched antibody pairs: MP01528-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01528-2, TMI GD2/CD28H Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84727-1-PBS. Detection antibody: 84727-3-PBS. Standard: Eg2264. Range: 0.781-25 ng/mL.