

For Research Use Only

NUAK2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01626-3

Capture Antibody Information

Catalog Number:
84849-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242220B3

Reactivity:
human

Immunogen Catalog Number:
Ag35302

Conjugate:
Unconjugated

Full name:
NUAK family, SNF1-like kinase, 2

Gene ID:
81788

Detection Antibody Information

Catalog Number:
84849-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242220C1

Reactivity:
human

GenBank:
NM_030952

Immunogen Catalog Number:
Ag35302

Conjugate:
Unconjugated

Full name:
NUAK family, SNF1-like kinase, 2

Gene ID:
81788

Applications

Tested Applications:
Sandwich ELISA

Range:
0.195-12.5 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01626-3 targets NUAK2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: NUAK2 Recombinant antibody, PBS Only (Capture) 84849-2-PBS (242220B3). 100 µg. Concentration 1 mg/ml.

Detection antibody: NUAK2 Recombinant antibody, PBS Only (Detector) 84849-4-PBS (242220C1). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01626-3, Human NUAk2 Recombinant Matched Antibody Pair - PBS only. 84849-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag35302. 84849-4-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL.

The mean NUAk2 concentration was determined to be 14.95 ng/mL in A549 cell extract based on a 2.3 mg/mL extract load.