

For Research Use Only

# PTH Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01717-3

## Capture Antibody Information

Catalog Number:  
84960-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242608F8

Reactivity:  
human

Conjugate:  
Unconjugated

Full name:  
parathyroid hormone

Gene ID:  
5741

## Detection Antibody Information

Catalog Number:  
84960-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242608F7

Reactivity:  
human

GenBank:  
NM\_000315.4

Conjugate:  
Unconjugated

Full name:  
parathyroid hormone

Gene ID:  
5741

## Applications

Tested Applications:  
Sandwich ELISA

Range:  
9.77-625 pg/mL (Sandwich ELISA)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01717-3 targets PTH in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: PTH Recombinant antibody, PBS Only (Capture) 84960-4-PBS (242608F8). 100 µg. Concentration 1 mg/ml.

Detection antibody: PTH Recombinant antibody, PBS Only (Detector) 84960-3-PBS (242608F7). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Sandwich ELISA standard curve of MP01717-3, Human PTH Recombinant Matched Antibody Pair - PBS only. 84960-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3191. 84960-3-PBS was HRP conjugated as the detection antibody. Range: 9.77-625 pg/mL



Serum of six individual healthy human donors was measured. The PTH concentration of detected samples was determined to be 44.3 pg/mL with a range of 28.9-54.3 pg/mL.



Plasma of six individual healthy human donors was measured. The PTH concentration of detected samples was determined to be 54.2 ng/mL with a range of 32.7-102.4 pg/mL.