

For Research Use Only

PDLIM1,CLP36 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01749-2

Capture Antibody Information

Catalog Number:
85027-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242585E4

Reactivity:
human

Immunogen Catalog Number:
Ag2274

Conjugate:
Unconjugated

Full name:
PDZ and LIM domain 1

Gene ID:
9124

Detection Antibody Information

Catalog Number:
85027-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242585C8

Reactivity:
human

GenBank:
BC000915

Immunogen Catalog Number:
Ag2274

Conjugate:
Unconjugated

Full name:
PDZ and LIM domain 1

Gene ID:
9124

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01749-2 targets PDLIM1,CLP36 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PDLIM1,CLP36 Recombinant antibody, PBS Only (Capture) 85027-2-PBS (242585E4). 100 µg.
Concentration 1 mg/ml.

Detection antibody: PDLIM1,CLP36 Recombinant antibody, PBS Only (Detector) 85027-1-PBS (242585C8). 100 µg.
Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01749-2, PDLIM1/CLP36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85027-2-PBS. Detection antibody: 85027-1-PBS. Standard: Ag2274. Range: 0.313-40 ng/mL.